

FIXED

Product/Brand	Composition details
Nonprecious	Argeloy NPV <i>Argen</i> © 72% Ni, 15% Cr, 9% Mo, 2% Al, 1.8% Be
NP Tilitte	Talladium <i>Tilitte</i> ® 76% Ni, 13.5% Cr, 6% Mo, 4% Ti
Noble	Argelite 58P <i>Argen</i> © 57.9% Pd, 30% Ag, 6% Sn, 4% In, 2% Zn
High Noble	Argedent Euro 40 <i>Argen</i> © 40% Au, 39.4% Pd, 10% Ag, 7.7% In, 1.4% Ga
Premium	Argedent 75 <i>Argen</i> © 75.1 Au, 12% Pd, 10% Ag, 1.8% In, 1.09% Sn
Captek™	Captek™ Nano 86% Au, 4.4% Pt, 6.3% Pd, 2.8% Ag, .5% Ir
Zirconia	TDS Fusion >99.5% ZrO2+ HfO2 , <5% Al2O3, <5% Other

FULL CAST

Product/Brand	Composition details
Nonprecious	Argeloy NPV <i>Argen</i> © 72% Ni, 15% Cr, 9% Mo, 2% Al, 1.8% Be
Noble silver	Argelite 58P <i>Argen</i> © 57.9% Pd, 30% Ag, 6% Sn, 4% In, 2% Zn
Noble gold	Argenco Y+ <i>Argen</i> © 2% Au, 34.9% Pd, 30% Ag, 30% In
High Noble	Argenco 40HN <i>Argen</i> © 40% Au, 20% Pd, 20% Ag, 17% In, 3% Zn

REMOVABLES

Product/Brand	Composition details
Framework I Nobileum	<i>NobilStar</i> ™ 65% Co, 27.7% Cr, 5.5% Mo
Framework II Nobileum	<i>NobilStar Ultra</i> ™ 64% Co, 28.5% Cr, 5.25% Mo

PORCELAINS

Brand	Usage
VMK® 95	Vita PFM, Captek, & Jacket Crowns with Vita Shades
VM® 9	Vita PFZirconia with Vita shades
IPS Inline®	Ivoclar PFM & Captek with Chromoscope shades
e.max®	Ivoclar Pressed ceramic crowns, veneers, & PFZ w/ Chromoscope Shades

STONES

Brand	Description
Die & Opposing	ETI Global Direct (California) Ivory, type 4
Base Stone	ETI Global Direct (California) Blue, flowable (poured with vibration)
Mounting Stone	ETI Global Direct (California) Plaster White



All materials have passed the United States FDA 501 (k) certification standards dating back to 1996 under Title 21, Chapter 1, Subchapter H of the United States Code of Federal Regulations.

Mabel Dental Lab assures you that the health and safety of our consumers comes first.

FDA Registration: 3009522026